

شهادة اعتماد

رقم ML 005



يقر نظام الاعتماد العراقي بأن:
مختبر النظائر المشعة
العراق- بغداد- الحارثية

تم اعتماده وفقا لمتطلبات المواصفة ISO 15189:2012
(المختبرات الطبية- متطلبات الكفاءة والجودة)

في مجال:
اختبارات الكيمياء السريرية

شرط التوافق مع متطلبات المواصفة اعلاه ومتطلبات IQAS الخاصة بالاعتماد
مجال الاعتماد المرفق بالشهادة يعتبر جزءا لا يتجزأ منها
يمكن الحصول على الاصدار الاحدث من مجال الاعتماد من خلال الموقع الالكتروني
<http://www.iraqaccreditation.net>

يكون الاعتماد نافذا من ٢٠٢٠/٠٤/٣٠ الى ٢٠٢٢/٠٤/٢٩
تاريخ منح الاعتماد لأول مرة
٢٠٢٠/٠٤/٣٠

د. نوري صباح الدليمي
وزير التخطيط

عبد الواحد محمد ابراهيم
معاون مدير عام الهيئة العراقية للاعتماد

Ministry of planning
Iraqi Organization for Accreditation
IQAS

ACCREDITATION CERTIFICATE

No.ML 005



Iraqi Accreditation System Certify that:

Al- Nadaer Clinical Laboratory

Iraq –Baghdad- Al-Harhiya

Is accredited according to the requirements of the standard ISO 15189:2012
(Medical laboratories — Requirements for quality and competence)

In the field of:

Clinical Chemistry Testing

This accreditation is subject to with the above standard & IQAS requirements
The scope of accreditation is attached to the certificate & considered as part of it
The most recent issue of the accreditation scope is available on IQAS website

<http://www.iraqaccreditation.net>

Accreditation is valid From 30/04/2020 To 29/04/2022

Initial accreditation date
30/04/2020

Abdulwahid M. Ibrahim
Deputy General Manager of
Iraqi Organization for Accreditation

Dr. Noori Sabah AL-Doulimi
Minister of Planning

	استمارة مجال الاعتماد Scope of Accreditation form	نظام الاعتماد العراقي IQAS
Organization address: Iraq –Baghdad- Al-Harthiya	Organization name : Al- Nadaer Clinical Laboratory	Accreditation no.: ML 005
Signature : 	Accreditation is valid From 30/04/2020 To 29/04/2022	Issue no. : 01

Branch/ Sector	Parameters/Test	Type of sample or system material	Principle of measurement	Instrument	Units
Clinical chemistry	Glucose	Serum	Hexokinase	Randox RX imola	mg/dl
	Cholesterol		Cholesterol Oxidase - Abell Kendall		mg/dl
	Triglyceride		Lipase/GPO-PAP no correction		mg/dl
	Creatinine		Enzymatic UV method (340nm)		mg/dl
	Urea		Urease, kinetic		mg/dl
	Calcium		Arsenazo		mg/dl
	Phosphate inorganic		Phosphomolybdate UV		mg/dl
	Alt		Tris buffer without P5P		U/l
	Ast		Tris buffer without P5P		U/l
	Uric Acid (urate)		Uricase Peroxidase with ascorb. ox @ 546nm		mg/dl